

MCU™

Multi-Cervical Unit

SITE REQUIREMENT SPECIFICATIONS



MCU

Height: 77 in / 196 cm

Width: 31 in / 79 cm

Depth: 47 in / 119 cm

Computer Cart w/Monitor

Height: 59 in / 150 cm

Width: 20 in / 51 cm

Depth: 30 in / 76 cm

Electrical Requirements

- ◆ Dedicated Circuit
- ◆ US/Canada: 110-120 Volts,
50-60 Hz, 14 Amps
- ◆ International: 220-240 Volts,
50-60 Hz, 7 Amps